

## Collecting a Drinking Water Monitoring Program CHEMISTRY Sample

### Sample Containers

Use **two** containers:

- One polyethylene screw-capped bottle capable of holding at least 250mL water, and
- One 250mL polyethylene screw-capped bottle with added acid preservative (5mL Nitric Acid 35% w/w with 0.14mL Hydrochloric Acid 17.5% w/w) as supplied by FASS.

### Caution: Acid Preservative for Metals in Water

The acid preservative is corrosive. It can cause severe burns and eye damage and may cause respiratory irritation.

A Safety Data Sheet (SDS) is supplied with the acidified bottles. Read this SDS before using the bottles.

A Safe Work Practice (SWP) for handling the acidified bottles and collecting the sample is attached to these sampling instructions. Read the SWP before sampling.

### Label


- Use a NSW Health barcoded label for the correct supply system, current year and sample type.

### Allocated Chemistry sample

- Use a **yellow** label with a barcode beginning with 7.
- Using an indelible pen record the **site code**, **time** and **date** of collection on the larger label.
- If the system is fluoridated conduct a fluoride test on a separate sample and record the **fluoride field result** on the label.
- Place the label evenly around the non-acidified sample bottle (not the lid) so that the entire barcode can be scanned at the laboratory.
- Attach the smaller matching barcoded label to the acidified bottle.

### Repeat and Additional Samples

- Repeat sample: use a **pink** label with a barcode starting with 6. Additional sample: use a **blue** label with a barcode starting with 8.
- Record the **site code**, **time and date** of collection and **fluoride field result** on the label.
- Place the label evenly around the non-acidified sample bottle.
- Record the supply system code (2 letters followed by 2 numbers), and the sampling date and time on a plain label and attach it to the acidified bottle.






## Collecting the Sample

- ◆ Flush the lines for at least 3 minutes before collecting the samples.
- ◆ Fill the non-acidified bottle to the brim.
- ◆ For the acidified bottle:
  - ◆ **Do not** rinse the bottle prior to taking the sample.
  - ◆ Fill the bottle to the brim, taking care **not to overfill**. Avoid splashing.

## Packaging the Samples as Dangerous Goods

- ◆ Place the non-reactive absorbent mat/material into an insulated cooler.
- ◆ Place the samples into the cooler with sufficient freezer bricks to keep the samples cool ( $\leq 6^{\circ}\text{C}$ ) during transport. Ensure the acidified samples are packed so that they remain upright in the cooler and are prevented from moving during transportation.
- ◆ Include a copy of the Safety Data Sheet in the cooler.
- ◆ Pack the cooler into a box strong enough to withstand conditions normally encountered during transport.
- ◆ Attach a copy of the FASS Address Label to the outside of the box.
- ◆ Attach a “Dangerous Goods” (orange) sticker, a “Corrosive” sticker and a “This Way Up” sticker to the outside of the box.
- ◆ Complete a Dangerous Goods consignment form.
- ◆ Complete a Sender’s Declaration form for Dangerous Goods. The following information must be included:
  - ◆ Shipping name - Corrosive liquid, Acidic, Inorganic, NOS
  - ◆ Technical Name – Mixture of nitric acid and hydrochloric acid
  - ◆ UN Number – 3264
  - ◆ DG Class – 8
  - ◆ Packing Group – II
  - ◆ Sender’s full name address and telephone number
  - ◆ FASS address:
 

480 Weeroona Rd,  
Lidcombe NSW 2141
- ◆ Dispatch the samples to FASS as soon as possible.

<b>URGENT WATER SAMPLES</b>		
Delivery Address		
Forensic & Analytical Science Service 480 Weeroona Road LIDCOMBE NSW 2141		
Laboratory	<input type="checkbox"/> Microbiology	<input type="checkbox"/> Chemistry / Organic Chemistry
After hours telephone number	0419 215 490	0413 984 105
<b>FRAGILE</b>		

